

Ravencourt  
**LIVING**

How to use your  
**Rosebud Reminder Clock**




Designed with health professionals especially  
for people living with memory problems

[www.ravencourtliving.com](http://www.ravencourtliving.com)

# Contents

What is the Rosebud Reminder Clock? .....	Pg3
Getting started.....	Pg3
Choosing the right clock display .....	Pg4
What to do when clocks change? .....	Pg4
The settings menu.....	Pg5-6
How to set reminder? .....	Pg7
What reminders are in my Rosebud?.....	Pg8
Personalising your Rosebud Reminder Clock.....	Pg9
Planning your reminders .....	Pg10
Customer feedback .....	Pg11



# What is the Rosebud Reminder Clock?

The Rosebud Reminder Clock was designed with health professionals specifically for people living with memory problems.

The Rosebud is a mains powered digital clock with

- 4 different clock displays
- 20 helpful reminders to aid with daily routine.

Rosebud Reminder is an image or text that appears on the digital screen along with a voice prompt.

Rosebud Reminder clock works well beside the bed, in the living room, kitchen and beside a calendar. You may have more than one Rosebud throughout the house.

## Getting started



Peel off the protected film from the screen.






Plug the power cable into the DC IN socket on the right side of the frame. Plug the cable into the mains and switch on at the socket.

Your Rosebud Reminder Clock is now ready to use

## Choosing the right calendar mode

When clock has been switched on you can choose one of 4 clock displays by pressing OK button. Each time you press OK button display will show one of calendar modes cycling through options.



## What to do when clocks change?

Your Rosebud Reminder clock will automatically move forwards 1 hour in spring and backwards by 1 hour in autumn in the line with UK daylight saving

# The Settings Menu

Set Time	< 07:10 >
Set Date	< 03-01-2019 >
Time Mode	< 24 Hours >
Date Mode	< Day-Month-Year >
Brightness - day	< 10 >
Brightness - night	< 1 >
Volume	< 6 >

## Set Time:

- Press OK button then LEFT or RIGHT button to underline the hour <07:10> change the hour by using UP or DOWN button.
- Press LEFT or RIGHT button to underline the minutes <07:10> change minutes by using UP or DOWN button.
- Press OK button to confirm the time

### **Rosebud times of day:**

- Morning: 5am to 12noon
- Afternoon: 12noon to 5pm
- Evening: 5pm to 10pm
- Night: 10pm to 5am

## Set Date:

- Press OK button then press LEFT or RIGHT button to underline the Day, Month or Year <03-01-2019> change the day, month or year by using UP or DOWN button.  
<03-01-2019> - change day  
<03-01-2019> - change month  
<03-01-2019> - change year
- Press OK button to confirm the date

## Time Mode:

- Press LEFT or RIGHT buttons to change between 24 hour or 12 hour time format, eg.23:00 or 11:00.

## Date Mode:

- Press LEFT or RIGHT button to choose one of two options <Day-Month-Year> OR <Month-Day-Year>

## Brightness – day / Brightness – night:

- Press LEFT or RIGHT button to change the brightness level for daytime and night time. The screen will change in brightness to help you make the right choice. The larger the number the brighter gets screen.


## Volume:

- Press LEFT or RIGHT button to change the volume from 1 to 10. The higher the number the louder the reminder.





# How to set a reminder?



- Press OK button twice to highlight the list of reminders
- Press UP and DOWN buttons to highlight the reminder and press OK button to confirm.
- Press the RIGHT button to underline the time and use UP and DOWN buttons to amend the time for your reminder.
- To choose the frequency of your reminder Press the RIGHT button again to underline Daily, Weekly, Future Date, Off or Once.
- For Weekly reminders press OK when Weekly is underlined and use UP or DOWN to choose the correct day. Press OK followed by MENU button to confirm.
- For Future date reminders press OK when Future date is underlined and use UP or DOWN to choose the date. Press OK followed by MENU button to confirm.

# What Reminders are already in my Rosebud?

**It's time for your lunch**



**It's time for your television programme**



**You have a hairdressers appointment today**



**Your carer is visiting today**

**It's time for a cup of tea**



**Your family are visiting today**

**Check the front door is locked**



**It's time for a cup of coffee**



**You have a doctors appointment today**



**It's time to feed your cat**



**Your home helper is visiting today**

**It's time for bed**



**It's time for your breakfast**



**It's time for your dinner**



**It's time to feed your dog**



**It's time to take your medication**



**It's time for a glass of water**



**Your lunch is in the fridge**



**Would you like to use the toilet?**



At your chosen time and day the built in reminders will first play a multi tonal chime. Then the image with text will appear on the display and the reminder will be spoken "You have a doctor appointment today". The reminder will repeat for 1 minute. You can stop the reminder at any time by pressing the OK button.



# Personalising your Rosebud Reminder Clock

1. First you must create your text or picture, the easiest way to do this is by using a program such as Microsoft Paint (you can use other similar programs). Select image size 1024 pixels wide X 768 pixels high. Then create your chosen image and save.
2. Create a folder on the chosen USB and name the folder, eg, Remember to feed the fish.
3. Drag and drop your chosen image into the folder you have created.
4. Record the audio message that you want in your reminder. This is easy to do on a smart phone or tablet and then emailed to your PC.
5. Drag and drop your sound file into the folder you have created. Without sound your reminder will not play. The length of the audio determines the length of the reminder. To repeat the audio, drop the audio file into the folder multiple times.
6. Each folder on the USB is a separate reminder.
7. For video reminders film a video in landscape and drop this into a folder on the USB card and name the folder appropriately.
8. Plug the USB card into the appropriate slot on your Rosebud reminder clock, you can then select from your own list of reminders as well as our pre-set list.



# Customer feedback

We are constantly striving to improve our products and the service that we offer. We also design new and innovative products for elderly care, memory loss and low vision. Any new product ideas or feedback on how to improve our products or services will be gratefully received. Please send your feedback and ideas to me directly for free by addressing an envelope with just the following “**Freepost RAVENCOURT**” if the envelope has any other address details on then it will not be delivered. I truly look forward to hearing from you.

Kind Regards

Oli Story

Director

Ravencourt Ltd



If you have any problems setting up your Rosebud or need advice on how to use your Rosebud, please contact us directly.

**info@ravencourtliving.com**

**01780-489100**

 **Ravencourt**  
**LIVING**

**The Rosebud Reminder Clock**

For people living with memory problems

**Wednesday afternoon**

**15:33**

**5 September 2018**

**It's now  
Thursday  
afternoon**



**Wednesday afternoon**

**Wednesday afternoon**

**4 different clock displays**

**20 pre-set reminders – or create your own**

Thank you for your purchase

**[info@ravencourtliving.com](mailto:info@ravencourtliving.com)**

**01780-489100**